REF **5050501** 5x1 ml
 LOT **1646CH** 04/2010
 IVD For in vitro medical device
 Per uso diagnostico in vitro

Uso

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Bibliografia

Department of Labor, Occupational Safety and Health Administration 29 CFR Part 1910.1030 Fed. Register July 1, 1998;6:267-280
 Council Directive 90/679/EEC. Official journal of the European directives N°L374 from Dec.31,1990:1-12

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Use

This product is intended for in quality control of Direct HDL-C, Direct LDL-C, Lipoprotein (a), Apolipoprotein A-1, Apolipoprotein B, Apolipoprotein A-II, Apolipoprotein C-II, Apolipoprotein C-III, Apolipoprotein E, cholesterol and triglycerides .

Summary

Lyophilized control serum prepared from human serum.
 The concentrations and activities of the components have been adjusted for a wide range of automatic analyzers and the assigned values are the means of at least 30 replicate determinations using various clinical chemistry systems.

Preparation

Open one bottle of control avoiding the loss of lyophilizate.
 Reconstitute the lyophilized serum with exactly 1.0 ml of bi-distilled water.
 Close the bottle and let it stand for 30 minutes out of bright light.
 Before to use mix by inverting the vial.
 Do not shake the vial as the formation of foam should avoid.

Storage and stability

Prior to reconstitution, when stored at 2-8 °C, the control is stable until the expire date.

After the reconstitution the serum is stable:

7 days	if stored at	2-8°C
4 weeks	if frozen at	-20° C

Caution

APO A-II APO C-III	28 days	if stored at	2-8°C
APO C-II	14 days	if stored at	2-8°C
APO E	1 day	if stored at	2-8°C
LP (a)	16 weeks	if frozen at	-20° C

Values may drop by up to 10 % for Direct LDL Cholesterol when stored for 4 weeks at -20° C.

Note

- Avoid direct light exposure, contamination and evaporation.
- In case of compliant or quality control request, refer to the lot number on the package or to the lot number of the single vial.
- Non properly reconstitution and / or storage can feed erratic results.

Precaution in use

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998).

Human serum is obtained using only blood of donors tested by FDA method and found non-reactive for HbsAg and negative for antibodies to HIV-1/2 and HCV.
 However, no known test method can offer complete assurance that product derived from Human blood transmit infectious diseases, this product should be handled as potentially infectious biological material with caution according to good laboratory practice.

Waste management

Please refer to the local legal requirements.

Bibliography

Department of Labour, Occupational Safety and Health Administration 29 CFR Part 1910.1030 Fed. Register July 1, 1998;6:267-280
 Council Directive 90/679/EEC. Official journal of the European directives N°L374 from Dec.31,1990:1-12

Symbols

	Biological risk
	CE Mark (98/79 CE regulation)
	in vitro medical device
	Batch Code
	Use by
	Storage temperature limits
	Read instruction for use
	Manufacturer in UK

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