REF 705CN 6x1 ml
IVD For in vitro medical device
 Per uso diagnostico in vitro

Uso

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Use

Normal control, for in vitro use, for monitoring the quality control of results obtained in the measuring of the antistreptolysin-0, C- reactive protein and rheumatoid factors in turbidimetry. Product from human origin.

Calibration values

The concentrations and activities of the components have been adjusted for a wide range of automatic analyzers and the assigned values are specific for each lot and they are listed below:

| BATH 200509 | MEAN | RANGE |
|--------------------|------|----------|
| ASO UI/ml | 113 | 79 – 147 |
| CRP mg/l | 13 | 9 – 17 |
| RF UI/ml | 30 | 21 – 39 |

Preparation

Ready to use

Storage and stability

When stored at 2-8 °C, it is stable until the expire date.

Note

- In case of compliant or quality control request, refer to the lot number on the package or to the lot number of the single vial.
- Not properly storage and / or contamination can feed not correct results.

Precaution in use

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998).

 Product from human origin is obtained using only blood of donors tested by FDA method and found non-reactive for HbsAg and negative for antibodies to HIV-1/2 and HCV. However, no known test method can offer complete assurance that product derived from Human blood transmit infectious diseases, this product should be handled as potentially infectious biological material with caution according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contact with skin, eyes and mucous membranes.



Waste management

Please refer to the local legal requirements.

Bibliography

Biosafety in Microbiological and Biomedical Laboratories. CDC/NIH manual, 4th edition, 1999.

Symbols

| | |
|---|-------------------------------|
|  | Biological risk |
|  | CE Mark (98/79 CE regulation) |
|  | in vitro medical device |
|  | Batch Code |
|  | Use by |
|  | Storage temperature limits |
|  | Read instruction for use |
|  | Manufacturer in Spain |

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